



Updated: September 4, 2020

Downloaded from www.dabir-surfaces.com/cleaners

Approved Cleaners (or equivalent):

Product Name	Manufacturer
CaviWipes	Metrex
Clorox Healthcare® Bleach Germicidal Wipes	The Clorox Company
Clorox® Dispatch Bleach Germicidal Cleaner	The Clorox Company
Expose® II 256	Diversey, Inc.
MICRO-KILL™ Bleach Germicidal Wipes	Medline Industries, Inc.
MICRO-KILL™+ Disinfectant Wipes	Medline Industries, Inc.
Oxivir® 1 Wipes	Diversey, Inc.
Oxivir® Five 16 Concentrate (US) (properly diluted)	Diversey, Inc.
Oxivir® TB Wipes Based	Diversey, Inc.
Oxycide™ Daily Disinfectant Cleaner	Ecolab Inc.
PeridoxRTU® Sporicidal Disinfectant and Cleaner	Contec, Inc.
Rubbing Alcohol (99%)	Not applicable
Sani-Cloth® Bleach Germicidal Disposable Wipes	Nice-Pak/PDI, Inc.
Sani-Cloth® Plus Germicidal Disposable Cloth	Nice-Pak/PDI, Inc.
Super Sani-Cloth® Germicidal Wipe	Nice-Pak/PDI, Inc.
Virex® II 256 One-Step Disinfectant Cleaner and Deodorant	Diversey, Inc.

Unapproved Cleaners*:

Product Name	Manufacturer
Clorox Commercial Solutions® Clorox® Disinfecting Wipes	Clorox Professional Products Company
Lysol® Brand II Disinfecting Wipes	Reckitt Benckiser Inc.
MICRO-KILL® One Germicidal Alcohol Wipes	Medline Industries, Inc.
Quat Plus TB	Medline Industries, Inc.
Sani-Cloth® AF3 Germicidal Disposable Wipes	Nice-Pak/PDI, Inc.
Virex® TB	Diversey, Inc.

*Test results have shown that the above "unapproved cleaners" will reduce Dabir® product life with repeated use.

NOTE: For more information on cleaning / disinfecting Dabir® products, visit:

<http://www.dabir-surfaces.com/IFU>